

SMALL VOLUME GLASSWARE

ROSETT GLASS COOLING CELLS

The Rosett cooling cell enables uniform treatment at low temperatures. The cell is placed in an ice bath. The ultrasonic energy forces the sample to circulate repeatedly under the probe and throughout the cooling arms.

30 ml Rosett cooling cell.

Part No. 830-00003

300 ml Rosett cooling cell.

Part No. 830-00001



GLASS COOLING CELLS

10 ml cooling cell with water jacket.

Part No. 830-00009

100 ml cooling cell with water jacket.

Part No. 830-00010

